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## Reliable screening test for novel Influenza A (H1N1)

**SD BIOLINE Influenza Ag rapid test can detect 2009 new strain of Influenza A(H1N1) virus isolated from human patients in Korea.**

#### Tested By

D.H Kwon, DVM  
Division of Influenza and Respiratory Viruses  
National Institute of Health  
Korea Centers for Disease Control and Prevention

#### Tested viruses and Concentration

2 cases of 2009 Influenza A /H1N1 viruses isolated from human patients in Korea.

- Sample No. 1 : 1.78X10<sup>4</sup> TCID<sub>50</sub>/ml
- Sample No. 2 : 4.57X10<sup>5</sup> TCID<sub>50</sub>/ml

#### Test Method

Each cultured virus was 10-fold diluted and 1 : 1 mixed with 100 µl of lysis buffer in test tube. The strip was dipped into the test tube and the results was interpreted within 15-20 min.

#### Test Results (reported on May 22, 2009)

virus	Concentration of diluted virus							
	10 <sup>-1</sup>	10 <sup>-2</sup>	10 <sup>-3</sup>	10 <sup>-4</sup>	10 <sup>-5</sup>	10 <sup>-6</sup>	10 <sup>-7</sup>	10 <sup>-8</sup>
No. 1	+	+	(weak)	-	-	-	-	-
No. 2	+	+	(weak)	+/-	-	-	-	-

### Certificate of Analytical Sensitivity

Product Name : SD BIOLINE Influenza Ag (Standard Diagnostics, Inc.)

We, BioNote, Inc. hereby certify that [SD BIOLINE Influenza Ag] can detect all A type influenza virus as follows. All tests were carried out at BioNote, and collaborating institutes in the world.

JinSik Oh  
Ph.D., DVM  
Executive Director  
BioNote, Inc

H type	Subtype	Strain	Disease
H1	H1N1	A/Duck/Totoni/723/80	Avian Influenza
	H1N1	A/PR/8/34	Avian Influenza
	H1N1	A/Swine/Korea/GC0503/2005	Swine Influenza
	H1N2	A/Swine/Korea/GC0502/2005	Swine Influenza
H2	H2N2	A/Duck/Hokkaido/17/01	Avian Influenza
	H2N2	A/Singapore/1/57	Avian Influenza
H3	H3N2	A/Swine/Korea/GC0407/2005	Swine Influenza
	H3N2	A/Swine/Korea/Suwon/107/1999	Swine Influenza
	H3N2	A/Swine/TX98	Swine Influenza
	H3N2	A/Canine/Korea/01/07	Canine Influenza
	H3N2	A/Canine/Korea/02/07	Canine Influenza
	H3N2	A/Canine/Korea/03/07	Canine Influenza
	H3N2	A/Canine/China/?	Canine Influenza

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#### Clinical Evaluation 1

Source: J Lab Med Qual Assur 2005;27:243-249  
(Full text article enclosed)

#### Evaluation of Two Lateral-Flow Chromatographic Membrane Immunoassays for Rapid Detection of Influenza Virus in Limited Respiratory Specimens

Le Thi Quynh Mai\*, Pham Thi Hien\*, Nguyen Le Khanh Hang \*, J S Oh \*\*, G W Ha\*\*, J A Kwon\*\*\*, C K Lee\*\*\*, K N Lee\*\*\*

\* Department of Virology, National Institute of Hygiene and Epidemiology, No1 Yersin Str, Hanoi 10,000 Vietnam  
\*\* Anigen Biotechnology Institute, Woncheon-dong, Yeonatong-gu, Suwon-si, Korea  
\*\*\*Department of Laboratory Medicine, Guro Hospital, Korea University, Seoul, Korea

**Background:** The diagnosis of influenza based on clinical grounds alone may be inaccurate, because the presenting symptoms of influenza are similar to those caused by other infectious agents. We evaluate two influenza rapid tests, SD BIOLINE Influenza Ag (Standard Diagnostic Inc., Yongin, Korea) and QuickVue™ Influenza Test (Quidel corporation, San Diego, USA) with influenza virus culture and RT-PCR.

**Methods:** The two commercially available rapid test kits, SD BIOLINE Influenza Ag and QuickVue™ Influenza Test, for influenza virus detection were evaluated with 189 respiratory specimens collected during Dec. 2004 to Nov. 2005 in Vietnam and compared with viral culture and RT-PCR.

**Results:** Overall, the SD BIOLINE Influenza Ag and QuickVue™ Influenza Test showed high sensitivities (88.4% and 82.6%, respectively) and high specificities (99.0% and 99.0%, respectively), high positive predictive value (PPV) (98.7% and 98.6%, respectively) and high negative predictive value (NPV) (91.1% and 87.2%, respectively).

**Conclusion:** Both SD BIOLINE Influenza Ag and QuickVue™ Influenza Test were easy to perform and showed high sensitivity and can be used as an additional tool for rapid diagnosis of influenza virus.

**Key words:** Rapid test, Influenza virus, Sensitivity, Specificity

Table. Viruses used for the analytical reactivity

Viral strain	Viral type(origin)	Viral strain	Viral type(origin)
A/PR/8/34	A (Human)	A1/Denver/1/57	A (Human)
A1/FM/1/47	A (Human)	Alice (H3N2)A	A
A/Mal/302/54 (H1N1)	A (Human)	MRC-2	A (Human)
A/WSN/33 (H1N1)	A (Human)	B/Hong Kong/5/72	B (Human)
A/Port Chalmers/1/73	A	B/Lee/40	B (Human)
A/Victoria/3/75	A (Human)	B/Taiwan/2/62	B (Human)
A/Weiss/43	A (Human)	B/Russia/69	B (Human)
* A/Swine/1976/31 (H1N1)	A (Hog)	B/GL/1739/54	B (Human)
H2S-15	A	B/R75	B (Human)
A2/Aichi/2/68 (H3N2)	A (Human)	B/Mass/3/66	B (Human)
A/Hong Kong/8/68 (H3N2)	A (Human)	B/Allen/45	B (Human)
A/NWS/33	A (Human)	B/Maryland/1/59	B (Human)
A/Equine 1/Prague/1/56	A (Horse)	B/rit	B (Human)
A/Equine 2/Miami/1/63	A (Horse)	B/R5	B (Human)
* A/Swine/Iowa/15/90	A (Hog)	B/R22 Barbara	B (Human)

### Results

#### Reactivity test of rapid test kits in reference panel

All 18 influenza virus A and 12 influenza B virus gave positive reactions for both two rapid tests. All 60 other microorganisms, including 25 viruses and 35 bacteria other than influenza virus type A or B tested for cross-reactivity were negative by both two rapid tests.

#### Clinical Evaluation 2

Source: Cvi, Aug. 2007, p. 1050-1052 Vol. 14, No. 8  
(Full text article enclosed)

CLINICAL AND VACCINE IMMUNOLOGY, Aug. 2007, p. 1050-1052  
1556-8811/07/8(8):00+0 doi:10.1128/CVI.00465-06  
Copyright © 2007, American Society for Microbiology. All Rights Reserved.

Vol. 14, No. 8

#### Clinical Evaluation of the SD Bioline Influenza Virus Antigen Test for Rapid Detection of Influenza Viruses A and B in Children and Adults during the Influenza Season<sup>V</sup>

Young Yoo,<sup>1</sup> Jang Wook Sohn,<sup>2</sup> Dae Won Park,<sup>2</sup> Jeong Yeon Kim,<sup>2</sup> Hye Kyung Shin,<sup>4</sup> Yoon Lee,<sup>1</sup> Ji Tae Choung,<sup>1</sup> Chang Kyu Lee,<sup>3</sup> and Min Ja Kim<sup>2,4\*</sup>

<sup>1</sup>Department of Pediatrics, <sup>2</sup>Division of Infectious Diseases, Department of Internal Medicine, <sup>3</sup>and Department of Laboratory Medicine, <sup>4</sup>College of Medicine, and Institute of Emerging Infectious Diseases, Korea University, Seoul, Korea

Received 6 December 2006/Returned for modification 13 March 2007/Accepted 4 June 2007

The performance of the SD Bioline rapid antigen test kit for influenza virus detection was evaluated with 295 respiratory specimens during the influenza season. The overall sensitivity and specificity of the SD Bioline test were 61.9% and 96.8% for the influenza A virus antigen and 54.5% and 100% for the influenza B virus antigen, respectively. The results were consistent with peak influenza activities.

Table. Performance of two rapid tests for the detection of influenza viruses

Rapid test and virus detected	No. of specimens that were:				Sensitivity (%)	Specificity (%)	Predictive value (%)	
	Rapid test positive, culture positive	Rapid test positive, culture negative	Rapid test negative, culture positive	Rapid test negative, culture negative			Positive	Negative
<b>SD Bioline test</b>								
Influenza A virus	26	8	16	245	61.9	96.8	76.5	93.1
Influenza B virus	18	0	15	262	54.5	100	94.6	94.9
QuickVue test	35	12	40	208	46.7	94.5	74.5	83.9

### Clinical Evaluation 3

(Full text article enclosed)

#### Evaluation of 19 Commercial Rapid Antigen (Point-of-Care) Assays in Detecting Influenza A H5N1, H3N2, H1N1 and B viruses.

Taylor J. McPhie K, Druce J, Birch C, Dwyer DE\*

Centre for Infectious Diseases & Microbiology Laboratory Services, Institute of Clinical Pathology and Medical Research, Westmead Hospital, Westmead NSW, and Victorian Infectious Diseases Reference Laboratory, North Carlton, Victoria, Australia

Table 1. Analytical sensitivity depending on various influenza virus type

	SD BIOLINE Influenza Antigen
H5N1(y irradiated) Indon E5	10 <sup>-2</sup>
H5N1(y irradiated) HN3028	10 <sup>-2</sup>
H3N2 (y irradiated) A/Fujian/114/2002-like	10 <sup>-2</sup>
H3N2 A/Fujian/114/2002-like	10 <sup>-2</sup>
H1N1(y irradiated) A/New Caledonia/20/99-like	10 <sup>-2</sup>
H1N1 A/New Caledonia/20/99-like	10 <sup>-2</sup>
B B/Shanghai/361/2002	10 <sup>-2</sup>
Volume of virus used	100 µl

### Clinical Evaluation 4

(Full text article enclosed)

ADDRESS : 227-5, KUGAL-DONG GIHEUNG-GU  
YONGIN-SI KYUNGGI-DO KOREA  
TEL : 82.31.283.3423 FAX : 82.31.281.8814  
E-mail : system@gcwp.co.kr

Project title: Analytical sensitivity of SD BIOLINE Influenza Ag (Standard Diagnostics, Inc.)

Evaluator: DaeSub Song, DVM, PhD.

Green Cross Veterinary Products Ltd.

Approved by *Kjongman*  
Head of Research Institute of  
Green Cross Veterinary Products Ltd.

Title *Head of Research Institute  
of GCVP*

Table 1. HA titers of swine influenza viruses

virus	HA titer
A/Swine/Korea/GC050312005(H 1 N 1)	32 HAU
A/Swine/Korea/GC0502/2005(H 1 N2)	8 HAU
3A/Swine/Korea/GC0407/2005(H3N2)	16 HAU

#### Results

No false positive was observed. The kit could detect 3 swine influenza viruses HI N 1, H 1 N2, and H3N2. Sensitivity-wise, the kit appeared to detect at least 1132 HA unit in a sample when tested using a SIVs isolate.

#### Conclusion

The kit appears to be capable of detecting swine influenza viruses

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
### Now it is used in Mexico to screen new Influenza A(H1N1)

#### Effort to stop Influenza outbreaks

Since new Influenza A H1N1 caused epidemic outbreaks, the whole world have put all efforts into stopping further transmission of new virus. As a screening test, SD BIOLINE Influenza antigen has been used all over the world. We've offered 1.3 million Influenza Ag rapid tests for Mexico through official bid program. Now it's your turn to use our rapid test.

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**CONTROL TÉCNICO Y REPRESENTACIONES, S.A. DE C.V.**
SERVILOO, S.A.